

adding a latex reagent directly to the sample +

- 7 to reacting the hemolysed whole blood sample in an agglutination reaction to
- 8 form a reaction product wherein a predetermined antigen in the hemolysed whole blood
- 9 sample specifically reacts with an antibody immobilized onto an insoluble carrier to ^{of the latex reagent}
- 10 provide the reaction product;
- 11 irradiating the reaction product in the sample with radiation which
- 12 includes a wavelength range which is substantially free from absorption by both
- 13 hemoglobin and the hemolysis reagent; and
- 14 measuring only in the wavelength range which is substantially free from
- 15 absorption by both hemoglobin and the hemolysis reagent, an absorbance of the incident
- 16 radiation by the reaction product to determine the quantity of antigens in the sample.
-

Q1
cont.

Q1
cont.